

For Research Use Only

# SPHK2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68871-2-PBS

## Basic Information

Catalog Number: 68871-2-PBS	GenBank Accession Number: BC010671	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 56848	CloneNo.: 2D11E4
Source: Mouse	UNIPROT ID: Q9NRA0	
Isotype: IgG1	Full Name: sphingosine kinase 2	
Immunogen Catalog Number: AG33745	Calculated MW: 618 aa, 65 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50264-1, SPHK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68871-1-PBS. Detection antibody: 68871-2-PBS. Standard:null. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50264-2, SPHK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68872-1-PBS. Detection antibody: 68871-2-PBS. Standard:Ag10569. Range: 0.195-25 ng/mL.