

For Research Use Only

SPHK2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68872-1-PBS

Basic Information

Catalog Number: 68872-1-PBS	GenBank Accession Number: BC010671	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 56848	CloneNo.: 3D12B9
Source: Mouse	UNIPROT ID: Q9NRA0	
Isotype: IgG1	Full Name: sphingosine kinase 2	
Immunogen Catalog Number: AG10569	Calculated MW: 618 aa, 65 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50264-2, SPHK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68872-1-PBS. Detection antibody: 68871-2-PBS. Standard: Ag10569. Range: 0.195-25 ng/mL