

For Research Use Only

IGFBP7 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68874-4-PBS

Basic Information

Catalog Number:

68874-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_001553.2

GeneID (NCBI):

3490

UNIPROT ID:

Q16270-1

Full Name:

insulin-like growth factor binding protein 7

Calculated MW:

29 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C4D3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50270-2, IGF1BP7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68874-3-PBS. Detection antibody: 68874-4-PBS. Standard: Eg0815. Range: 3.125-200 ng/mL.