

For Research Use Only

BRAF Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68875-1-PBS

Basic Information

Catalog Number:

68875-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28574

GenBank Accession Number:

BC101757

GeneID (NCBI):

673

UNIPROT ID:

P15056

Full Name:

v-raf murine sarcoma viral oncogene

homolog B1

Calculated MW:

766 aa, 84 kDa

Purification Method:

Protein A purification

CloneNo.:

1F5E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50272-1, BRAF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68875-1-PBS. Detection antibody: 68875-2-PBS. Standard:Ag28574. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50272-2, BRAF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68875-1-PBS. Detection antibody: 68875-3-PBS. Standard:Ag28574. Range: 0.391-100 ng/mL.