

For Research Use Only

APRIL, TNFSF13 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68876-2-PBS

Basic Information

Catalog Number:

68876-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG26281

GenBank Accession Number:

BC008042

GeneID (NCBI):

8741

UNIPROT ID:

O75888

Full Name:

tumor necrosis factor (ligand)
superfamily, member 13

Calculated MW:

27 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E10B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50273-1, APRIL/TNFSF13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68876-1-PBS. Detection antibody: 68876-2-PBS. Standard: Ag26281. Range: 1.563-100 ng/mL.