

For Research Use Only

# APRIL, TNFSF13 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68876-2-PBS



## Basic Information

Catalog Number:

68876-2-PBS

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG26281

GenBank Accession Number:

BC008042

GeneID (NCBI):

8741

UNIPROT ID:

O75888

Full Name:

tumor necrosis factor (ligand)  
superfamily, member 13

Calculated MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

2E10B5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50273-1, APRIL/TNFSF13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68876-1-PBS. Detection antibody: 68876-2-PBS. Standard: Ag26281. Range: 1.563-100 ng/mL.