

For Research Use Only

# APP Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68880-1-PBS

## Basic Information

Catalog Number:

68880-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22408

GenBank Accession Number:

BC065529

GeneID (NCBI):

351

UNIPROT ID:

P05067

Full Name:

amyloid beta (A4) precursor protein

Purification Method:

Protein G purification

CloneNo.:

5G5H3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50278-1, APP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68880-1-PBS. Detection antibody: 68880-2-PBS. Standard:Ag22408. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50278-2, APP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68880-1-PBS. Detection antibody: 68880-3-PBS. Standard:Ag22408. Range: 6.25-100 ng/mL.