

For Research Use Only

APP Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68880-2-PBS

Basic Information

Catalog Number: 68880-2-PBS	GenBank Accession Number: BC065529	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 351	CloneNo.: 3H1G8
Source: Mouse	UNIPROT ID: P05067	
Isotype: IgG1	Full Name: amyloid beta (A4) precursor protein	
Immunogen Catalog Number: AG22408		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50278-1, APP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68880-1-PBS. Detection antibody: 68880-2-PBS. Standard: Ag22408. Range: 3.125-100 ng/mL.