

For Research Use Only

CENPJ Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68881-1-PBS

Basic Information

Catalog Number: 68881-1-PBS	GenBank Accession Number: BC024209	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 55835	CloneNo.: 2F4F11
Source: Mouse	UNIPROT ID: Q9HC77	
Isotype: IgG1	Full Name: centromere protein J	
Immunogen Catalog Number: AG28390	Calculated MW: 1338 aa, 153 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50279-1, CENPJ Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68881-1-PBS. Detection antibody: 68881-2-PBS. Standard:Ag28390. Range: 6.25-200 ng/mL