

For Research Use Only

DLL4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68882-3-PBS

Basic Information

Catalog Number: 68882-3-PBS	GenBank Accession Number: BC106950	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 54567	CloneNo.: 1F3D12
Source: Mouse	UNIPROT ID: Q9NR61	
Isotype: IgG2a	Full Name: delta-like 4 (Drosophila)	
Immunogen Catalog Number: AG17483	Calculated MW: 685 aa, 75 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50281-2, DLL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68882-1-PBS. Detection antibody: 68882-3-PBS. Standard: Ag16140. Range: 0.391-12.5 ng/mL