

For Research Use Only

BCL2A1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68883-3-PBS

Basic Information

Catalog Number: 68883-3-PBS	GenBank Accession Number: BC016281	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 597	CloneNo.: 2F6D11
Source: Mouse	UNIPROT ID: Q16548	
Isotype: IgG2b	Full Name: BCL2-related protein A1	
Immunogen Catalog Number: AG34180	Calculated MW: 175 aa, 20 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50284-2, BCL2A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68883-3-PBS. Detection antibody: 68883-1-PBS. Standard:Ag34180. Range: 0.195-200 ng/mL.