

For Research Use Only

BCL2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68887-1-PBS

Basic Information

Catalog Number:

68887-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG27874

GenBank Accession Number:

BC027258

GeneID (NCBI):

596

UNIPROT ID:

P10415

Full Name:

B-cell CLL/lymphoma 2

Calculated MW:

26 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2F3G9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50290-1, BCL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68887-1-PBS. Detection antibody: 60178-2-PBS. Standard: Ag3508. Range: 1.563-200 ng/mL.