

For Research Use Only

Integrin alpha-E Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68892-1-PBS

Basic Information

Catalog Number:

68892-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG17892

GenBank Accession Number:

BC113436

GeneID (NCBI):

3682

UNIPROT ID:

P38570

Full Name:

integrin, alpha E (antigen CD103,
human mucosal lymphocyte antigen
1; alpha polypeptide)

Calculated MW:

1179 aa, 130 kDa

Purification Method:

Protein A purification

CloneNo.:

1E1C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50294-1, Integrin alpha-E Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68892-1-PBS. Detection antibody: 68892-2-PBS. Standard:null. Range: null.



Cytometric bead array standard curve of MP50294-1, Integrin alpha-E Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68892-1-PBS. Detection antibody: 68892-2-PBS. Standard:Ag17892. Range: 3.125-100 ng/mL.