

For Research Use Only

# Integrin alpha-E Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68892-2-PBS

## Basic Information

Catalog Number:

68892-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17892

GenBank Accession Number:

BC113436

GeneID (NCBI):

3682

UNIPROT ID:

P38570

Full Name:

integrin, alpha E (antigen CD103,  
human mucosal lymphocyte antigen  
1; alpha polypeptide)

Calculated MW:

1179 aa, 130 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C10G1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50294-1, Integrin alpha-E Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68892-1-PBS. Detection antibody: 68892-2-PBS. Standard:null. Range: null.



Cytometric bead array standard curve of MP50294-1, Integrin alpha-E Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68892-1-PBS. Detection antibody: 68892-2-PBS. Standard:Ag17892. Range: 3.125-100 ng/mL.