

For Research Use Only

ACTA2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68895-2-PBS

Basic Information

Catalog Number: 68895-2-PBS	GenBank Accession Number: BC017554	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 59	CloneNo.: 1B12C9
Source: Mouse	UNIPROT ID: P62736	
Isotype: IgG1	Full Name: actin, alpha 2, smooth muscle, aorta	
Immunogen Catalog Number: AG19478	Calculated MW: 377 aa, 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50296-2, ACTA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68895-2-PBS. Detection antibody: 67735-2-PBS. Standard:Ag6088. Range: 3.125-50 ng/mL