

For Research Use Only

FOXA1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68896-3-PBS

Basic Information

Catalog Number: 68896-3-PBS	GenBank Accession Number: BC033890	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3169	CloneNo.: 1D9F1
Source: Mouse	UNIPROT ID: P55317	
Isotype: IgG2a	Full Name: forkhead box A1	
Immunogen Catalog Number: AG14502	Calculated MW: 473 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50297-2, FOXA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68896-1-PBS. Detection antibody: 68896-3-PBS. Standard: Ag14502. Range: 1.563-100 ng/mL.