

For Research Use Only

ICT1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68898-2-PBS

Basic Information

Catalog Number:

68898-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33109

GenBank Accession Number:

BC015335

GeneID (NCBI):

3396

UNIPROT ID:

Q14197

Full Name:

immature colon carcinoma transcript

1

Calculated MW:

24 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C6E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50300-1, ICT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68898-1-PBS. Detection antibody: 68898-2-PBS. Standard: Ag33109. Range: 3.125-200 ng/mL.