

For Research Use Only

GTF2IRD1 Monoclonal antibody

Catalog Number: 68899-1-Ig



Basic Information

Catalog Number: 68899-1-Ig	GenBank Accession Number: BC018136	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 9569	CloneNo.: 1E5H3
Source: Mouse	UNIPROT ID: Q9UHL9	Recommended Dilutions: WB 1:2000-1:10000
Isotype: IgG2a	Full Name: GTF2I repeat domain containing 1	
Immunogen Catalog Number: AG10761	Calculated MW: 959 aa, 106 kDa	
	Observed MW: 125 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : SH-SY5Y cells, hTERT-RPE1 cells, HeLa cells, TT cells, MKN-45 cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68899-1-Ig (GTF2IRD1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.