

For Research Use Only

TRIM11 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68904-1-PBS

Basic Information

| | | |
|-------------------------------------|--|---|
| Catalog Number: 68904-1-PBS | GenBank Accession Number: BC011629 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 81559 | CloneNo.: 2D6D1 |
| Source: Mouse | UNIPROT ID: Q96F44 | |
| Isotype: IgG1 | Full Name: tripartite motif-containing 11 | |
| Immunogen Catalog Number: AG1308 | Calculated MW: 52 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50306-1, TRIM11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68904-1-PBS. Detection antibody: 68904-2-PBS. Standard:Ag1308. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50306-2, TRIM11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68904-1-PBS. Detection antibody: 68904-3-PBS. Standard:Ag1308. Range: 0.781-100 ng/mL