

For Research Use Only

PMS2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68905-2-PBS

Basic Information

Catalog Number:

68905-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG34237

GenBank Accession Number:

BC093921

GeneID (NCBI):

5395

UNIPROT ID:

P54278

Full Name:

PMS2 postmeiotic segregation
increased 2 (*S. cerevisiae*)

Calculated MW:

862 aa, 96 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3B11C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50308-1, PMS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68905-1-PBS. Detection antibody: 68905-2-PBS. Standard: Ag34237. Range: 1.563-200 ng/mL.