

For Research Use Only

PMS2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68905-3-PBS

Basic Information

Catalog Number: 68905-3-PBS	GenBank Accession Number: BC093921	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 5395	CloneNo.: 3E7E7
Source: Mouse	UNIPROT ID: P54278	
Isotype: IgG1	Full Name: PMS2 postmeiotic segregation increased 2 (S. cerevisiae)	
Immunogen Catalog Number: AG34237	Calculated MW: 862 aa, 96 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50308-2, PMS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68905-3-PBS. Detection antibody: 68905-4-PBS. Standard: Ag34237. Range: 1.563-100 ng/mL.