

For Research Use Only

HYPK Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68906-4-PBS



Basic Information

Catalog Number: 68906-4-PBS	GenBank Accession Number: BC019262	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 25764	CloneNo.: 1H1C3
Source: Mouse	UNIPROT ID: Q9NX55	
Isotype: IgG2b	Full Name: Huntingtin interacting protein K	
Immunogen Catalog Number: AG34843	Calculated MW: 129 aa, 15 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50312-3, HYPK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68906-4-PBS. Detection antibody: 68906-5-PBS. Standard: Ag34843. Range: 0.391-100 ng/mL.