

For Research Use Only

SENP1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68907-2-PBS

Basic Information

Catalog Number: 68907-2-PBS	GenBank Accession Number: BC045639	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 29843	CloneNo.: 3C7E1
Source: Mouse	UNIPROT ID: Q9POU3	
Isotype: IgG2a	Full Name: SUMO1/sentrin specific peptidase 1	
Immunogen Catalog Number: AG18708	Calculated MW: 644 aa, 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50313-1, SENP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68907-1-PBS. Detection antibody: 68907-2-PBS. Standard: Ag18708. Range: 0.391-100 ng/mL.