

For Research Use Only

# LAMA4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68908-1-PBS

## Basic Information

|                                      |                                       |   |
|--------------------------------------|---------------------------------------|---|
| Catalog Number:<br>68908-1-PBS       | GenBank Accession Number:<br>BC004241 | Purification Method:<br>Protein A Magarose purification |
| Size:<br>1 mg/ml                     | GeneID (NCBI):<br>3910                | CloneNo.:<br>3F3D11                                     |
| Source:<br>Mouse                     | UNIPROT ID:<br>Q16363                 |   |
| Isotype:<br>IgG2b                    | Full Name:<br>laminin, alpha 4        |   |
| Immunogen Catalog Number:<br>AG35773 | Calculated MW:<br>13 kDa              |   |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50314-1, LAMA4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68908-1-PBS. Detection antibody: 68908-2-PBS. Standard: Ag35773. Range: 1.563-100 ng/mL.