

For Research Use Only

ROPN1L Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68909-1-PBS

Basic Information

Catalog Number: 68909-1-PBS	GenBank Accession Number: BC014607	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 83853	CloneNo.: 1C11D3
Source: Mouse	UNIPROT ID: Q96C74	
Isotype: IgG1	Full Name: ropporin 1-like	
Immunogen Catalog Number: AG34015	Calculated MW: 230 aa, 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50320-1, ROPN1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68909-1-PBS. Detection antibody: 68909-2-PBS. Standard: Ag34015. Range: 0.781-100 ng/mL.