

For Research Use Only

# ROPN1L Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68909-2-PBS

## Basic Information

Catalog Number: 68909-2-PBS	GenBank Accession Number: BC014607	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 83853	CloneNo.: 1F3B9
Source: Mouse	UNIPROT ID: Q96C74	
Isotype: IgG1	Full Name: ropporin 1-like	
Immunogen Catalog Number: AG34015	Calculated MW: 230 aa, 26 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50320-1, ROPN1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68909-1-PBS. Detection antibody: 68909-2-PBS. Standard: Ag34015. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50320-2, ROPN1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68909-3-PBS. Detection antibody: 68909-2-PBS. Standard: Ag34015. Range: 0.781-100 ng/mL.