

For Research Use Only

CCDC46 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68910-4-PBS

Basic Information

Catalog Number: 68910-4-PBS	GenBank Accession Number: BC131805	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 201134	CloneNo.: 2F12F5
Source: Mouse	UNIPROT ID: Q8N8E3	
Isotype: IgG1	Full Name: coiled-coil domain containing 46	
Immunogen Catalog Number: AG34261	Calculated MW: 955 aa, 113 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50322-2, CCDC46 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68910-3-PBS. Detection antibody: 68910-4-PBS. Standard:Ag34261. Range: 0.781-100 ng/mL.