

For Research Use Only

TSLP Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68911-2-PBS

Basic Information

Catalog Number: 68911-2-PBS	GenBank Accession Number: BC016720	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 85480	CloneNo.: 4G8B11
Source: Mouse	UNIPROT ID: Q969D9	
Isotype: IgG2a	Full Name: thymic stromal Lymphopoietin	
Immunogen Catalog Number: AG4771	Calculated MW: 15 kDa, 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50325-1, TSLP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68911-1-PBS. Detection antibody: 68911-2-PBS. Standard: Ag4771. Range: 0.098-100 ng/mL.