

For Research Use Only

WT1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68912-2-PBS

Basic Information

Catalog Number: 68912-2-PBS	GenBank Accession Number: BC032861	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 7490	CloneNo.: 3B8H1
Source: Mouse	UNIPROT ID: P19544	
Isotype: IgG1	Full Name: Wilms tumor 1	
Immunogen Catalog Number: AG19541	Calculated MW: 449 aa, 49 kDa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50330-1, WT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68912-1-PBS. Detection antibody: 68912-2-PBS. Standard: Ag19541. Range: 6.25-200 ng/mL