

For Research Use Only

# CD28 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68913-1-PBS

## Basic Information

<b>Catalog Number:</b> 68913-1-PBS	<b>GenBank Accession Number:</b> BC093698	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 940	<b>CloneNo.:</b> 1H2C5
<b>Source:</b> Mouse	<b>ENSEMBL Gene ID:</b> ENSG00000178562	
<b>Isotype:</b> IgG1	<b>UNIPROT ID:</b> P10747	
<b>Immunogen Catalog Number:</b> AG24058	<b>Full Name:</b> CD28 molecule	
	<b>Calculated MW:</b> 220 aa, 25 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50331-1, CD28 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68913-1-PBS. Detection antibody: 68913-2-PBS. Standard: Ag24058. Range: 0.391-12.5 ng/mL