

For Research Use Only

CD28 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68913-2-PBS

Basic Information

Catalog Number:

68913-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG24058

GenBank Accession Number:

BC093698

GeneID (NCBI):

940

ENSEMBL Gene ID:

ENSG00000178562

UNIPROT ID:

P10747

Full Name:

CD28 molecule

Calculated MW:

220 aa, 25 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

2D1D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50331-1, CD28 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68913-1-PBS. Detection antibody: 68913-2-PBS. Standard:Ag24058. Range: 0.391-12.5 ng/mL