

For Research Use Only

CD28 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68913-3-PBS

Basic Information

Catalog Number: 68913-3-PBS	GenBank Accession Number: BC093698	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 940	CloneNo.: 1E12C5
Source: Mouse	ENSEMBL Gene ID: ENSG00000178562	
Isotype: IgG1	UNIPROT ID: P10747	
Immunogen Catalog Number: AG24058	Full Name: CD28 molecule	
	Calculated MW: 220 aa, 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50331-2, CD28 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68913-3-PBS. Detection antibody: 68913-4-PBS. Standard:Ag24058. Range: 1.563-200 ng/mL.