

For Research Use Only

# PARK7/DJ-1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68915-5-PBS

## Basic Information

Catalog Number:

68915-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28526

GenBank Accession Number:

BC008188

GeneID (NCBI):

11315

UNIPROT ID:

Q99497

Full Name:

Parkinson disease (autosomal  
recessive, early onset) 7

Calculated MW:

189 aa, 20 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4G9H9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50337-3, PARK7/DJ-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68915-4-PBS. Detection antibody: 68915-5-PBS. Standard: Ag28526. Range: 1.563-100 ng/mL.