

For Research Use Only

SLC4A2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68916-1-PBS



Basic Information

Catalog Number:

68916-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG23673

GenBank Accession Number:

BC028601

GeneID (NCBI):

6522

UNIPROT ID:

P04920

Full Name:

solute carrier family 4, anion
exchanger, member 2 (erythrocyte
membrane protein band 3-like 1)

Calculated MW:

1241 aa, 137 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B6A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50500-1, SLC4A2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68916-1-PBS. Detection antibody: 68916-2-PBS. Standard: Ag23673. Range: 0.781-100 ng/mL.