

For Research Use Only

KIF2A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68917-3-PBS

Basic Information

Catalog Number: 68917-3-PBS	GenBank Accession Number: BC031828	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3796	CloneNo.: 3C2A10
Source: Mouse	UNIPROT ID: O00139	
Isotype: IgG2b	Full Name: kinesin heavy chain member 2A	
Immunogen Catalog Number: AG34099	Calculated MW: 679 aa, 77 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50342-2, KIF2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68917-1-PBS. Detection antibody: 68917-3-PBS. Standard: Ag34099. Range: 0.195-100 ng/mL.