

For Research Use Only

PTCH1 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 68918-4-PBS

Basic Information

Catalog Number: 68918-4-PBS	GenBank Accession Number: NM_000264	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5727	CloneNo.: 3C12D5
Source: Mouse	UNIPROT ID: Q13635	
Isotype: IgG1	Full Name: patched homolog 1 (Drosophila)	
Immunogen Catalog Number: AG28106	Calculated MW: 161 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50345-2, PTCH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68918-3-PBS. Detection antibody: 68918-4-PBS. Standard: Ag28106. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50345-3, PTCH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68918-4-PBS. Detection antibody: 68918-5-PBS. Standard: Ag28106. Range: 0.391-100 ng/mL.