

For Research Use Only

DNMT3B Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68919-2-PBS

Basic Information

Catalog Number:

68919-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33907

GenBank Accession Number:

NM_001207055

GeneID (NCBI):

1789

UNIPROT ID:

Q9UBC3

Full Name:

DNA (cytosine-5-)-methyltransferase

3 beta

Calculated MW:

95 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A3C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50346-1, DNMT3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68919-1-PBS. Detection antibody: 68919-2-PBS. Standard: Ag33907. Range: 0.391-200 ng/mL.