

For Research Use Only

DNMT3B Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68919-4-PBS

Basic Information

Catalog Number: 68919-4-PBS	GenBank Accession Number: NM_001207055	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1789	CloneNo.: 2F10E2
Source: Mouse	UNIPROT ID: Q9UBC3	
Isotype: IgG1	Full Name: DNA (cytosine-5-)-methyltransferase	
Immunogen Catalog Number: AG33907	3 beta	
	Calculated MW: 95 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50346-2, DNMT3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68919-3-PBS. Detection antibody: 68919-4-PBS. Standard: Ag33907. Range: 0.391-200 ng/mL.