

For Research Use Only

# DNMT3B Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68919-6-PBS

## Basic Information

Catalog Number:

68919-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG33907

GenBank Accession Number:

NM\_001207055

GeneID (NCBI):

1789

UNIPROT ID:

Q9UBC3

Full Name:

DNA (cytosine-5-)-methyltransferase

3 beta

Calculated MW:

95 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F12B3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50346-4, DNMT3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68919-6-PBS. Detection antibody: 68919-5-PBS. Standard: Ag33907. Range: 0.391-25 ng/mL