

For Research Use Only

TMEM106B Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68920-1-PBS

Basic Information

| | | |
|---|---|---|
| Catalog Number: 68920-1-PBS | GenBank Accession Number: BC033901 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 54664 | CloneNo.: 1D12A3 |
| Source: Mouse | UNIPROT ID: Q9NUM4 | |
| Isotype: IgG2b | Full Name: transmembrane protein 106B | |
| Immunogen Catalog Number: AG14264 | Calculated MW: 31 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50326-4, TMEM106B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68920-1-PBS. Detection antibody: 68920-2-PBS. Standard: Ag14264. Range: 6.25-100 ng/mL.