

For Research Use Only

CCDC109A Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68926-3-PBS

Basic Information

Catalog Number: 68926-3-PBS	GenBank Accession Number: BC034235	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 90550	CloneNo.: 2G11E10
Source: Mouse	UNIPROT ID: Q8NE86	
Isotype: IgG1	Full Name: coiled-coil domain containing 109A	
Immunogen Catalog Number: AG22748	Calculated MW: 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50349-2, CCDC109A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68926-3-PBS. Detection antibody: 68926-4-PBS. Standard: Ag22748. Range: 3.125-100 ng/mL.