

For Research Use Only

NPC2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68927-3-PBS

Basic Information

Catalog Number: 68927-3-PBS	GenBank Accession Number: BC002532	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10577	CloneNo.: 2C7E3
Source: Mouse	UNIPROT ID: P61916	
Isotype: IgG2b	Full Name: Niemann-Pick disease, type C2	
Immunogen Catalog Number: AG34710	Calculated MW: 151 aa, 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50350-2, NPC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68927-3-PBS. Detection antibody: 68927-4-PBS. Standard: Ag34710. Range: 1.563-100 ng/mL.