

For Research Use Only

NPC2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68927-4-PBS

Basic Information

Catalog Number:

68927-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34710

GenBank Accession Number:

BC002532

GeneID (NCBI):

10577

UNIPROT ID:

P61916

Full Name:

Niemann-Pick disease, type C2

Calculated MW:

151 aa, 17 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D8D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50350-2, NPC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68927-3-PBS. Detection antibody: 68927-4-PBS. Standard: Ag34710. Range: 1.563-100 ng/mL.