

For Research Use Only

# IL28RA Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68930-1-PBS

## Basic Information

Catalog Number:

68930-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18830

GenBank Accession Number:

BC140872

GeneID (NCBI):

163702

UNIPROT ID:

Q8IU57

Full Name:

interleukin 28 receptor, alpha  
(interferon, lambda receptor)

Calculated MW:

520 aa, 58 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A12A3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50353-1, IL28RA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68930-1-PBS. Detection antibody: 68930-2-PBS. Standard: Ag18830. Range: 12.5-200 ng/mL