

For Research Use Only

IL28RA Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68930-2-PBS

Basic Information

Catalog Number: 68930-2-PBS	GenBank Accession Number: BC140872	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 163702	CloneNo.: 1D2F4
Source: Mouse	UNIPROT ID: Q8IU57	
Isotype: IgG3	Full Name: interleukin 28 receptor, alpha (interferon, lambda receptor)	
Immunogen Catalog Number: AG18830	Calculated MW: 520 aa, 58 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50353-1, IL28RA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68930-1-PBS. Detection antibody: 68930-2-PBS. Standard: Ag18830. Range: 12.5-200 ng/mL