

For Research Use Only

CD94 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68932-1-PBS

Basic Information

Catalog Number:

68932-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG24477

GenBank Accession Number:

BC028009

GeneID (NCBI):

3824

UNIPROT ID:

Q13241

Full Name:

killer cell lectin-like receptor
subfamily D, member 1

Calculated MW:

179 aa, 20 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3F4H10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50356-1, KLRD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68932-1-PBS. Detection antibody: 68932-2-PBS. Standard: Ag24477. Range: 0.781-100 ng/mL.