

For Research Use Only

# KLRD1 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 68932-2-PBS

## Basic Information

Catalog Number:

68932-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG24477

GenBank Accession Number:

BC028009

GeneID (NCBI):

3824

UNIPROT ID:

Q13241

Full Name:

killer cell lectin-like receptor  
subfamily D, member 1

Calculated MW:

179 aa, 20 kDa

Purification Method:

Protein G purification

CloneNo.:

2B10G12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50356-1, KLRD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68932-1-PBS. Detection antibody: 68932-2-PBS. Standard:Ag24477. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50356-2, KLRD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68932-2-PBS. Detection antibody: 68932-3-PBS. Standard:Ag24477. Range: 0.781-100 ng/mL.