

For Research Use Only

# LMAN1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68933-1-PBS

## Basic Information

Catalog Number: 68933-1-PBS	GenBank Accession Number: BC032330	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3998	CloneNo.: 3G2D10
Source: Mouse	UNIPROT ID: P49257	
Isotype: IgG1	Full Name: lectin, mannose-binding, 1	
Immunogen Catalog Number: AG4352	Calculated MW: 510 aa, 54 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50360-1, LMAN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68933-1-PBS. Detection antibody: 68933-2-PBS. Standard: Ag4352. Range: 3.125-100 ng/mL.