

For Research Use Only

PRKCQ Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68939-2-PBS

Basic Information

Catalog Number:

68939-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26331

GenBank Accession Number:

BC101465

GeneID (NCBI):

5588

UNIPROT ID:

Q04759

Full Name:

protein kinase C, theta

Calculated MW:

706 aa, 82 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D1G3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50373-1, PRKCC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68939-1-PBS. Detection antibody: 68939-2-PBS. Standard: Ag26331. Range: 0.781-25 ng/mL