

For Research Use Only

Integrin beta-6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68942-2-PBS



Basic Information

Catalog Number:

68942-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29015

GenBank Accession Number:

BC121178

GeneID (NCBI):

3694

UNIPROT ID:

P18564

Full Name:

integrin, beta 6

Calculated MW:

788 aa, 86 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3F4F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50377-1, Integrin beta-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68942-1-PBS. Detection antibody: 68942-2-PBS. Standard: Ag29015. Range: 0.391-50 ng/mL